

Constant-temperature water bath shaker

Features:

Widely applicable for laboratory researches on bacteria cultivation, fermentation, hybridization, chemical and biochemical reaction, enzymes and tissue research, which have a high requirement on precision of shaking speed and temperature. Timing Range: 1~9999min.

Specifications:

Model	DKZ-2	DKZ-3
Temperature Range	RT+5~99°C	
Temperature Control	0.1°C	
Temperature Uniformly	±0.1°C	
Shaking Frequency	30~150r/Min	
Shaking Speed Range	30~40mm Optional	
Input power	1250W	1650W
Interior Dimension (mm) W*D*H	438*310*250	618*310*250
Overall Dimension (mm) W*D*H	643*350*353	823*350*355

