

MCKD-05/MKCD-07 Mini Centrifuge



MCKD-05

- * with progressive high voltage technology,
- * speed won't be affected by voltage fluctuation
- * widely used and accept filter tubes
- * Complete with 2 types rotor and tube, optional with glass slide rotor.
- * low noise
- * double door design, more safety
- * measure up to global safety standard
- * with safety switch,
- * operation stop while opening the door
- * with rotor disassembly.

Model	MCKD-05	MCKD-07
MAX RPM	5000 RPM	7000 RPM
Max RCF	1500xg	2000xg
Rotor capacity	Rotor capacity:0.5m/1.5MI/2mL*8, 0.2mL PCR tubes*16 (double rows)	
Voltage	110V- 240V 50-60Hz	
Dimension(mm)	160x170x122mm	
Weight	0.5kg	