

SB-3200DTDN Ultrasonic Cleaner



Features:

- Large-screen LCD display
- Display the parameters setting menu
- Adjustable power:40%–100%
- Digital setting ultrasonic cleaning time
- Working time countdown display
- Display actual temperature of the cleaning tank
- Temperature setting Range: Ambient temperature–80°C
- Timing Range: 1-999 minutes adjustable
- Equipped with a special stainless steel basket, reduce sound cover
- The shell and cover are made of stainless steel
- The equipment controlled by Single-chip controller
- The last working parameters can be stored in memory after power-off

Main Technical Parameters:

Model	Inner Tank	Volume	Frequency	power	Heating power	Drainage	Rack
	L/W/H (mm)	L	KHz	W	W		
SB-120DTDP	240×140×100	3	40	120	150	/	√
SB-3200DTDN	300×155×150	6	40	180	150	/	√
SB-5200DTDN	300×240×150	10	40	200	800	√	√
SB-5200DTDN	300×240×150	10	40	250	800	√	√
SB-5200DTDN	300×240×150	10	40	300	800	√	√
SB25-12DTDN	500×300×150	22.5	40	500	1000	√	√
SB25-12DTDN	500×300×150	22.5	40	600	1000	√	√