

SB25-12YDTD Ultrasonic Cleaner

Main Features :

- LCD Background light screen
- Menu preference on screen
- Display working frequency
- Power of 40%-100% adjustable
- Digital display and setting ultrasonic cleaning time
- Display working countdown time
- Display actual temperature of the cleaning tank
- Temperature setting range: Ambient Temperature-80°C
- Timing range: 1-240 minutes adjustable
- Equipped with a special stainless steel basket, reduce sound cover
- Shell and cover are made of stainless steel
- Equipment controlled by Single-chip controller
- Last working parameters can be stored in memory after power-off



Main Technical Parameters:

Model	Inner Tank	Volume	Ultrasonic frequency	Ultrasonic power	Heating power	Drainage
	L/W/H (mm)	L	KHz	W	W	
SB-3200YDTD	300×155×150	6	40	180	120	/
SB-5200YDTD	300×240×150	10	40	300	320	√
SB-4200YDTD	300×240×200	15	40	400	320	√
SB25-12YDTD	500×300×150	22.5	40	500	600	√
SB25-12YDTD	500×300×200	30	40	600	800	√